

## PAR-4 (1-6) amide (human)

## Chemical Properties

CAS No. : 245443-51-0


Formula: C<sub>28</sub>H<sub>42</sub>N<sub>8</sub>O<sub>8</sub>

Molecular Weight: 618.68

Keep away from moisture

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.



## Biological Description

Description	PAR-4 (1-6) amide human, an N-terminal fragment of protease-activated receptor 4 (PAR4), induces platelet aggregation [1].
Targets(IC50)	Protease-activated Receptor
In vitro	PAR-4 (1-6) amide human (50, 100 μM) induces platelet aggregation [1].

## Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	1.6163 mL	8.0817 mL	16.1634 mL
5 mM	0.3233 mL	1.6163 mL	3.2327 mL
10 mM	0.1616 mL	0.8082 mL	1.6163 mL
50 mM	0.0323 mL	0.1616 mL	0.3233 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286

E\_mail: info@targetmol.com

Address: 34 Washington Street, Wellesley Hills, MA 02481