

PBX-7011

## Chemical Properties


CAS No. : 3003834-30-5

Formula: C<sub>24</sub>H<sub>21</sub>N<sub>3</sub>O<sub>6</sub>

Molecular Weight: 447.44

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.



## Biological Description

Description	PBX-7011, an active camptothecin derivative, binds to the DDX5 protein in cells and induces cell death through DDX5 protein degradation [1].
Targets(IC50)	Others, Survivin

## Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	2.2349 mL	11.1747 mL	22.3494 mL
5 mM	0.447 mL	2.2349 mL	4.4699 mL
10 mM	0.2235 mL	1.1175 mL	2.2349 mL
50 mM	0.0447 mL	0.2235 mL	0.447 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286 E\_mail: info@targetmol.com Address: 34 Washington Street, Wellesley Hills, MA 02481