


PI3K α -IN-16

Chemical Properties

CAS No. :	1448711-61-2
Formula:	C ₂₇ H ₂₄ N ₂ O ₂
Molecular Weight:	408.49
Storage:	Powder: -20°C for 3 years In solvent: -80°C for 1 year Actual storage temperature shall be subject to the COA.



Biological Description

Description	PI3K α -IN-16 (Z86) is a recently developed inhibitor of PI3K α , displaying an IC ₅₀ of 4.28 μ M. With potent efficacy, it inhibits the growth and migration of CRCs mediated by PI3K. Additionally, PI3K α -IN-16 also suppresses the Wnt signaling pathway [1].
Targets(IC ₅₀)	PI3K,Wnt/beta-catenin

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	2.448 mL	12.2402 mL	24.4804 mL
5 mM	0.4896 mL	2.448 mL	4.8961 mL
10 mM	0.2448 mL	1.224 mL	2.448 mL
50 mM	0.049 mL	0.2448 mL	0.4896 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481