


## PLK1/BRD4-IN-2

## Chemical Properties

CAS No. :	2251709-89-2
Formula:	C <sub>27</sub> H <sub>36</sub> N <sub>6</sub> O <sub>3</sub>
Molecular Weight:	492.61
Storage:	Powder: -20°C for 3 years   In solvent: -80°C for 1 year Actual storage temperature shall be subject to the COA.



## Biological Description

Description	PLK1/BRD4-IN-2 (compound 15), an analog of BI-2536, simultaneously inhibits both Polo-like kinase 1 (PLK1 IC <sub>50</sub> = 40 nM) and BRD4 bromodomain (BRD4-BD1 IC <sub>50</sub> = 28 nM) as a dual inhibitor [1].
Targets(IC <sub>50</sub> )	Epigenetic Reader Domain

## Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	2.030 mL	10.150 mL	20.300 mL
5 mM	0.406 mL	2.030 mL	4.060 mL
10 mM	0.203 mL	1.015 mL	2.030 mL
50 mM	0.0406 mL	0.203 mL	0.406 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286 E\_mail: info@targetmol.com Address: 34 Washington Street, Wellesley Hills, MA 02481