


PROTAC SMARCA2 degrader-3

Chemical Properties

CAS No. :	3024266-69-8
Formula:	C37H42N8O4
Molecular Weight:	662.78
Storage:	Keep away from direct sunlight Powder: -20°C for 3 years In solvent: -80°C for 1 year <small>Actual storage temperature shall be subject to the COA.</small>



Biological Description

Description	PROTAC SMARCA2 degrader-3, a degrader targeting the SWI/SNF ATPase subunits specifically SMARCA2, exhibits anticancer properties (WO2023244764A1; Compound 153) [1].
Targets(IC50)	Epigenetic Reader Domain,PROTACs

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	1.5088 mL	7.544 mL	15.088 mL
5 mM	0.3018 mL	1.5088 mL	3.0176 mL
10 mM	0.1509 mL	0.7544 mL	1.5088 mL
50 mM	0.0302 mL	0.1509 mL	0.3018 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481