


Rac1-IN-4

Chemical Properties

CAS No. :	2924486-45-1
Formula:	C ₂₃ H ₁₆ ClN ₅ O ₂
Molecular Weight:	429.86
Storage:	Powder: -20°C for 3 years In solvent: -80°C for 1 year Actual storage temperature shall be subject to the COA.



Biological Description

Description	Rac1-IN-4 (CAS: 2924486-45-1) is a Rac1 inhibitor [1].
Targets(IC50)	Ras

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	2.3263 mL	11.6317 mL	23.2634 mL
5 mM	0.4653 mL	2.3263 mL	4.6527 mL
10 mM	0.2326 mL	1.1632 mL	2.3263 mL
50 mM	0.0465 mL	0.2326 mL	0.4653 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481