

Resorufin acetate

Chemical Properties


CAS No. : 1152-14-3

Formula: C₁₄H₉NO₄

Molecular Weight: 255.23

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.



Biological Description

Description	Resorufin acetate, a fluorescent dye, is used to selectively detect hydrazine [1].
Targets(IC50)	Others

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	3.918 mL	19.5902 mL	39.1803 mL
5 mM	0.7836 mL	3.918 mL	7.8361 mL
10 mM	0.3918 mL	1.959 mL	3.918 mL
50 mM	0.0784 mL	0.3918 mL	0.7836 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286

E_mail:info@targetmol.com

Address:34 Washington Street,Wellesley Hills,MA 02481