


Solvent red 195

Chemical Properties

CAS No. :	72968-71-9
Formula:	C ₂₃ H ₂₉ N ₇ O ₄ S
Molecular Weight:	499.59
Storage:	Powder: -20°C for 3 years In solvent: -80°C for 1 year Actual storage temperature shall be subject to the COA.



Biological Description

Description	Solvent Red 195 is a red dye used in dyeing applications of plastic materials, along with its related products, such as Solvent Orange 60.
Targets(IC50)	Others

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	2.0016 mL	10.0082 mL	20.0164 mL
5 mM	0.4003 mL	2.0016 mL	4.0033 mL
10 mM	0.2002 mL	1.0008 mL	2.0016 mL
50 mM	0.040 mL	0.2002 mL	0.4003 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286 E_mail: info@targetmol.com Address: 34 Washington Street, Wellesley Hills, MA 02481