


Solvent yellow 124

Chemical Properties

CAS No. : 34432-92-3
 Formula: C₂₂H₃₁N₃O₂
 Molecular Weight: 369.5
 Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year
 Actual storage temperature shall be subject to the COA.



Biological Description

Description	Solvent Yellow 124 is a yellow dye used in the dyeing of plastic materials, along with related products like Solvent Orange 60.
Targets(IC50)	Others

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	2.7064 mL	13.5318 mL	27.0636 mL
5 mM	0.5413 mL	2.7064 mL	5.4127 mL
10 mM	0.2706 mL	1.3532 mL	2.7064 mL
50 mM	0.0541 mL	0.2706 mL	0.5413 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286 E_mail: info@targetmol.com Address: 34 Washington Street, Wellesley Hills, MA 02481