


## Solvent Yellow 98

## Chemical Properties

CAS No. : 12671-74-8  
 Formula: C<sub>36</sub>H<sub>45</sub>NO<sub>2</sub>S  
 Molecular Weight: 555.81  
 Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year  
 Actual storage temperature shall be subject to the COA.



## Biological Description

Description	Solvent Yellow 98, a yellow dye, belongs to a product series including Solvent Orange 60, which is utilized in the coloring of plastic materials.
Targets(IC50)	Others

## Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	1.7992 mL	8.9959 mL	17.9918 mL
5 mM	0.3598 mL	1.7992 mL	3.5984 mL
10 mM	0.1799 mL	0.8996 mL	1.7992 mL
50 mM	0.036 mL	0.1799 mL	0.3598 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481