

SOMCL-668

Chemical Properties


CAS No. : 1422251-09-9

Formula: C17H19NO

Molecular Weight: 253.34

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.



Biological Description

Description	SOMCL-668, a selective and potent sigma-1 receptor allosteric modulator, demonstrates positive modulation in ameliorating social deficits and cognitive impairment triggered by the selective sigma-1 agonist PRE084. Additionally, SOMCL-668 exhibits anti-seizure properties, making it useful for research in psychotic illnesses [1].
Targets(IC50)	Sigma receptor

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	3.9473 mL	19.7363 mL	39.4726 mL
5 mM	0.7895 mL	3.9473 mL	7.8945 mL
10 mM	0.3947 mL	1.9736 mL	3.9473 mL
50 mM	0.0789 mL	0.3947 mL	0.7895 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286

E_mail:info@targetmol.com

Address:34 Washington Street,Wellesley Hills,MA 02481