

SZ-3

Chemical Properties


CAS No. : 2556443-73-1

Formula: C₂₉H₄₃N₃O₃

Molecular Weight:

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.



Biological Description

Description

SZ-3, a cofilin inhibitor, reduces cofilin expression and protects SH-SY5Y cells from thrombin-induced cell death. Additionally, SZ-3 decreases microglial migration and proliferation and reduces inflammatory cytokine release [1].

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286

E_mail:info@targetmol.com

Address:34 Washington Street,Wellesley Hills,MA 02481