


TBK1/IKKε-IN-3

Chemical Properties

CAS No. : 851814-28-3
 Formula: C₁₆H₁₁NO₃
 Molecular Weight: 265.26
 Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year
 Actual storage temperature shall be subject to the COA.



Biological Description

Description	TBK1/IKKε-IN-3, a potent IP6K inhibitor, exhibits IC ₅₀ values of 3, 8.65, and 3.72 μM for IP6K1, IP6K2, and IP6K3, respectively. It is applicable in obesity disease research. [1]
Targets(IC ₅₀)	IκB/IKK

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	3.7699 mL	18.8494 mL	37.6989 mL
5 mM	0.754 mL	3.7699 mL	7.5398 mL
10 mM	0.377 mL	1.8849 mL	3.7699 mL
50 mM	0.0754 mL	0.377 mL	0.754 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286 E_mail: info@targetmol.com Address: 34 Washington Street, Wellesley Hills, MA 02481