


Thyroid hormone receptor beta agonist-1

Chemical Properties

CAS No. :	2140902-80-1
Formula:	C35H45O9P
Molecular Weight:	640.70
Storage:	Powder: -20°C for 3 years In solvent: -80°C for 1 year Actual storage temperature shall be subject to the COA.



Biological Description

Description	Compound 3, also known as Thyroid hormone receptor beta agonist-1, functions as a thyroid beta-agonist [1].
Targets(IC50)	Thyroid hormone receptor(THR)

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	1.5608 mL	7.804 mL	15.6079 mL
5 mM	0.3122 mL	1.5608 mL	3.1216 mL
10 mM	0.1561 mL	0.7804 mL	1.5608 mL
50 mM	0.0312 mL	0.1561 mL	0.3122 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481