


## TRRAP-IN-1

## Chemical Properties

CAS No. : 857254-94-5  
 Formula: C<sub>20</sub>H<sub>20</sub>O<sub>4</sub>  
 Molecular Weight: 324.37  
 Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year  
 Actual storage temperature shall be subject to the COA.



## Biological Description

Description	TRRAP-IN-1 (compound 1) acts as a potent inhibitor of the MYC:TRRAP interaction, which is significant in anti-tumor research [1] [2].
-------------	---

## Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	3.0829 mL	15.4145 mL	30.829 mL
5 mM	0.6166 mL	3.0829 mL	6.1658 mL
10 mM	0.3083 mL	1.5414 mL	3.0829 mL
50 mM	0.0617 mL	0.3083 mL	0.6166 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481