

Upifitamab rilsodotin

Chemical Properties

CAS No. : 2254119-00-9

Formula:

Molecular Weight:

Storage: -20°C for 1 year
Actual storage temperature shall be subject to the COA.

Biological Description

Description	Upifitamab rilsodotin (XMT-1536), an antibody-drug conjugate, specifically targets NaPi2b and is utilized in cancer research, particularly in studies related to ovarian cancer [1] [2].
Targets(IC50)	Sodium-dependent phosphate transporter

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481