

YIGSR

Chemical Properties

CAS No. : 110590-64-2


Formula: C₂₆H₄₂N₈O₈

Molecular Weight: 594.66

Keep away from moisture

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.



Biological Description

Description	YIGSR is a peptide that inhibits tumor growth and metastasis of leukemic cells [1].
Targets(IC50)	NO Synthase

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	1.6816 mL	8.4082 mL	16.8163 mL
5 mM	0.3363 mL	1.6816 mL	3.3633 mL
10 mM	0.1682 mL	0.8408 mL	1.6816 mL
50 mM	0.0336 mL	0.1682 mL	0.3363 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286 E_mail: info@targetmol.com Address: 34 Washington Street, Wellesley Hills, MA 02481