

KY371

Chemical Properties


CAS No. :

Formula: C₁₈H₂₉N₅O

Molecular Weight: 339.43

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.



Biological Description

Description	KY371, a click chemistry reagent containing an alkyne group, is a broad-range glycosidase inhibitor that can be used in conjunction with proteomic analysis of targeted proteins to study glycosidases [1].
Targets(IC50)	ADC Linker

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	2.9461 mL	14.7306 mL	29.4612 mL
5 mM	0.5892 mL	2.9461 mL	5.8922 mL
10 mM	0.2946 mL	1.4731 mL	2.9461 mL
50 mM	0.0589 mL	0.2946 mL	0.5892 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286 E_mail: info@targetmol.com Address: 34 Washington Street, Wellesley Hills, MA 02481