


## Mal-amido-PEG9-Val-Ala-PAB-SG3200

## Chemical Properties

CAS No. : 1817788-75-2  
 Formula: C75H101N9O24  
 Molecular Weight: 1512.65  
 Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year  
 Actual storage temperature shall be subject to the COA.



## Biological Description

Description	Mal-amido-PEG9-Val-Ala-PAB-SG3200 is a cleavable ADC linker conjugate employed in the production of antibody-drug conjugates (ADCs) [1].
Targets(IC50)	ADC Linker

## Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	0.6611 mL	3.3055 mL	6.6109 mL
5 mM	0.1322 mL	0.6611 mL	1.3222 mL
10 mM	0.0661 mL	0.3305 mL	0.6611 mL
50 mM	0.0132 mL	0.0661 mL	0.1322 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481