

## Glucocorticoid receptor agonist-3 Ala-Ala-Mal

## Chemical Properties


CAS No. : 3014393-35-9

Formula: C49H55FN4O12

Molecular Weight: 910.98

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.



## Biological Description

Description	Glucocorticoid Receptor Agonist-3 Ala-Ala-Mal (Compound Preparation 8) is a conjugate of an anti-human TNF $\alpha$ antibody and a glucocorticoid receptor agonist (GC). It is utilized in researching autoimmune and inflammatory diseases [1].
Targets(IC50)	Glucocorticoid Receptor, Drug-Linker Conjugates for ADC, TNF

## Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	1.0977 mL	5.4886 mL	10.9772 mL
5 mM	0.2195 mL	1.0977 mL	2.1954 mL
10 mM	0.1098 mL	0.5489 mL	1.0977 mL
50 mM	0.022 mL	0.1098 mL	0.2195 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286

E\_mail: info@targetmol.com

Address: 34 Washington Street, Wellesley Hills, MA 02481