


INX-P

Chemical Properties

CAS No. : 2734878-02-3
 Formula: C43H52BrN3O11
 Molecular Weight: 866.79
 Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year
 Actual storage temperature shall be subject to the COA.



Biological Description

Description	INX-P is a glucocorticosteroid-targeting antibody-drug conjugate (ADC) [1].
Targets(IC50)	Glucocorticoid Receptor,ADC Linker

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	1.1537 mL	5.7684 mL	11.5368 mL
5 mM	0.2307 mL	1.1537 mL	2.3074 mL
10 mM	0.1154 mL	0.5768 mL	1.1537 mL
50 mM	0.0231 mL	0.1154 mL	0.2307 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481