


## C3-Amide-C4-NH2

## Chemical Properties

CAS No. : 65389-82-4  
 Formula: C<sub>8</sub>H<sub>18</sub>N<sub>2</sub>O  
 Molecular Weight: 158.24  
 Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year  
 Actual storage temperature shall be subject to the COA.



## Biological Description

Description	C3-Amide-C4-NH2 is an effective linker [1].
Targets(IC50)	ADC Linker

## Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	6.3195 mL	31.5976 mL	63.1951 mL
5 mM	1.2639 mL	6.3195 mL	12.639 mL
10 mM	0.632 mL	3.1598 mL	6.3195 mL
50 mM	0.1264 mL	0.632 mL	1.2639 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481