


N-Boc-N-methyl-D-Valinol

Chemical Properties

CAS No. :	2165479-27-4
Formula:	C ₁₁ H ₂₃ NO ₃
Molecular Weight:	217.31
Storage:	Powder: -20°C for 3 years In solvent: -80°C for 1 year Actual storage temperature shall be subject to the COA.



Biological Description

Description	N-Boc-N-methyl-D-Valinol, an ADC linker with a BOC protecting group.
Targets(IC50)	ADC Linker

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	4.6017 mL	23.0086 mL	46.0172 mL
5 mM	0.9203 mL	4.6017 mL	9.2034 mL
10 mM	0.4602 mL	2.3009 mL	4.6017 mL
50 mM	0.092 mL	0.4602 mL	0.9203 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286 E_mail: info@targetmol.com Address: 34 Washington Street, Wellesley Hills, MA 02481