

## Anti-HN1 Antibody

### Chemical Properties

CAS No. :

Formula:

Molecular Weight:

Storage:

-20°C for 1 year

Actual storage temperature shall be subject to the COA.

### Biological Description

Description	Anti-HN1 Antibody is an inhibitory antibody that can be used in life science related research.
Targets(IC50)	Others

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286

E\_mail:info@targetmol.com

Address:34 Washington Street,Wellesley Hills,MA 02481