


(±)-Pellotine hydrochloride**Chemical Properties**

CAS No. : 26126-01-2
 Formula: C₁₃H₁₉NO₃·HCl
 Molecular Weight: 273.76
 Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year
 Actual storage temperature shall be subject to the COA.

**Biological Description**

Description	(±)-Pellotine has been found in <i>L. diffusa</i> and <i>L. fricii</i> , as well as <i>L. williamsii</i> . It decreases locomotor activity and the amount of rapid eye movement (REM) sleep in mice.
-------------	--

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	3.6528 mL	18.2642 mL	36.5283 mL
5 mM	0.7306 mL	3.6528 mL	7.3057 mL
10 mM	0.3653 mL	1.8264 mL	3.6528 mL
50 mM	0.0731 mL	0.3653 mL	0.7306 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286 E_mail: info@targetmol.com Address: 34 Washington Street, Wellesley Hills, MA 02481