

Butaverine

Chemical Properties


CAS No. : 55837-14-4

Formula: C₁₈H₂₇N₂O₂

Molecular Weight: 289.41

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.



Biological Description

Description	Butaverine is an antispasmodic compound.
Targets(IC50)	Others

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	3.4553 mL	17.2765 mL	34.5531 mL
5 mM	0.6911 mL	3.4553 mL	6.9106 mL
10 mM	0.3455 mL	1.7277 mL	3.4553 mL
50 mM	0.0691 mL	0.3455 mL	0.6911 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481