


Kv7.2 modulator 1

Chemical Properties

CAS No. :	3034884-49-3
Formula:	C ₂₂ H ₂₂ FN ₃ O ₂
Molecular Weight:	379.43
Storage:	Powder: -20°C for 3 years In solvent: -80°C for 1 year Actual storage temperature shall be subject to the COA.



Biological Description

Description	Kv7.2 Modulator 1 (compound 10), a modulator of the Kv7.2 channel, is utilized in epilepsy research.
Targets(IC50)	Potassium Channel

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	2.6355 mL	13.1777 mL	26.3553 mL
5 mM	0.5271 mL	2.6355 mL	5.2711 mL
10 mM	0.2636 mL	1.3178 mL	2.6355 mL
50 mM	0.0527 mL	0.2636 mL	0.5271 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481