


WP836

Chemical Properties

CAS No. : 266000-09-3
 Formula: C₃₁H₃₃N₃O₁₁
 Molecular Weight: 595.59
 Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year
Actual storage temperature shall be subject to the COA.



Biological Description

Description	WP836 is an effective antineoplastic agent, demonstrating cytotoxic properties.
Targets(IC50)	Others,DNA Alkylator/Crosslinker

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	1.679 mL	8.395 mL	16.7901 mL
5 mM	0.3358 mL	1.679 mL	3.358 mL
10 mM	0.1679 mL	0.8395 mL	1.679 mL
50 mM	0.0336 mL	0.1679 mL	0.3358 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481