

E3 ligase Ligand 27

Chemical Properties

CAS No. : 3034204-25-3


Formula: C₂₈H₃₈N₄O₆S

Molecular Weight: 558.69

Keep away from direct sunlight

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.



Biological Description

Description	E3 Ligase Ligand 27 (part of 15) serves as a ligand for the E3 ligase, employed in the creation of PROTAC Bcl-2/Bcl-xL Degrader-1.
Targets(IC50)	Ligands for E3 Ligase

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	1.7899 mL	8.9495 mL	17.899 mL
5 mM	0.358 mL	1.7899 mL	3.5798 mL
10 mM	0.179 mL	0.895 mL	1.7899 mL
50 mM	0.0358 mL	0.179 mL	0.358 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286

E_mail: info@targetmol.com

Address: 34 Washington Street, Wellesley Hills, MA 02481