


## Proxicromil

## Chemical Properties

CAS No. :	60400-92-2
Formula:	C <sub>17</sub> H <sub>18</sub> O <sub>5</sub>
Molecular Weight:	302.32
Storage:	Powder: -20°C for 3 years   In solvent: -80°C for 1 year Actual storage temperature shall be subject to the COA.



## Biological Description

Description	Proxicromil, a derivative of chromone, is an orally effective anti-asthmatic agent that does not exhibit hepatotoxicity.
Targets(IC50)	Others

## Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	3.3078 mL	16.5388 mL	33.0775 mL
5 mM	0.6616 mL	3.3078 mL	6.6155 mL
10 mM	0.3308 mL	1.6539 mL	3.3078 mL
50 mM	0.0662 mL	0.3308 mL	0.6616 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481