


Pan-RAS-IN-3

Chemical Properties

CAS No. :	3034588-80-9
Formula:	C ₃₈ H ₃₉ F ₃ N ₈ O ₃
Molecular Weight:	712.76
Storage:	Powder: -20°C for 3 years In solvent: -80°C for 1 year Actual storage temperature shall be subject to the COA.



Biological Description

Description	Pan-RAS-IN-3, a pan-Ras inhibitor, is applicable in the research of melanoma and acute myeloid leukemia.
Targets(IC50)	Ras

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	1.403 mL	7.015 mL	14.030 mL
5 mM	0.2806 mL	1.403 mL	2.806 mL
10 mM	0.1403 mL	0.7015 mL	1.403 mL
50 mM	0.0281 mL	0.1403 mL	0.2806 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286 E_mail: info@targetmol.com Address: 34 Washington Street, Wellesley Hills, MA 02481