

CSF1R-IN-24

Chemical Properties


CAS No. : 3024322-57-1

Formula: C₂₆H₂₄N₆O₃

Molecular Weight: 468.51

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.



Biological Description

Description	CSF1R-IN-24 (Example 134), an orally active inhibitor of CSF1R, effectively reduces the survival of human iPSC-derived microglia (hiPSC-MG).
Targets(IC50)	c-Fms

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	2.1344 mL	10.6721 mL	21.3443 mL
5 mM	0.4269 mL	2.1344 mL	4.2689 mL
10 mM	0.2134 mL	1.0672 mL	2.1344 mL
50 mM	0.0427 mL	0.2134 mL	0.4269 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286

E_mail: info@targetmol.com

Address: 34 Washington Street, Wellesley Hills, MA 02481