


S-palm P0(180-199)

Chemical Properties

CAS No. :	1240687-76-6
Formula:	C112H191N33O29S2
Molecular Weight:	2528.05
Storage:	Keep away from moisture Powder: -20°C for 3 years In solvent: -80°C for 1 year <small>Actual storage temperature shall be subject to the COA.</small>



Biological Description

Description	S-palm P0(180-199) is a peptide that can enhance MHC II class-restricted responses and is applicable in developing models of chronic inflammatory demyelinating polyneuropathy (CIDP).
-------------	--

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	0.3956 mL	1.9778 mL	3.9556 mL
5 mM	0.0791 mL	0.3956 mL	0.7911 mL
10 mM	0.0396 mL	0.1978 mL	0.3956 mL
50 mM	0.0079 mL	0.0396 mL	0.0791 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481