


HDAC2-IN-2

Chemical Properties

CAS No. : 332169-78-5
Formula: C₁₈H₁₅N₃O₃S
Molecular Weight: 353.40
Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year
Actual storage temperature shall be subject to the COA.



Biological Description

Description	HDAC2-IN-2 (compound 124) acts as an HDAC2 inhibitor, exhibiting a K _d value ranging from 0.1-1 μM.
Targets(IC50)	HDAC

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	2.8297 mL	14.1483 mL	28.2965 mL
5 mM	0.5659 mL	2.8297 mL	5.6593 mL
10 mM	0.283 mL	1.4148 mL	2.8297 mL
50 mM	0.0566 mL	0.283 mL	0.5659 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481