

## Genkwanol A

## Chemical Properties


CAS No. : 111103-90-3

Formula: C<sub>30</sub>H<sub>22</sub>O<sub>10</sub>

Molecular Weight:

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.



## Biological Description

Description

Genkwanol A was isolated from the roots of Wikstroemia indica.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286

E\_mail:info@targetmol.com

Address:34 Washington Street,Wellesley Hills,MA 02481