


## 3-Bromoisonicotinic acid

## Chemical Properties

CAS No. :	13959-02-9
Formula:	C <sub>6</sub> H <sub>4</sub> BrNO <sub>2</sub>
Molecular Weight:	202.007
Storage:	Powder: -20°C for 3 years   In solvent: -80°C for 1 year <small>Actual storage temperature shall be subject to the COA.</small>



## Biological Description

Description	3-Bromoisonicotinic acid, as the most important intermediate among isonicotinic acid derivatives, is mainly used to prepare anti tuberculosis drugs, as well as to synthesize amide, hydrazide, ester and other derivatives.
-------------	--

## Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	4.9502 mL	24.7512 mL	49.5025 mL
5 mM	0.990 mL	4.9502 mL	9.9005 mL
10 mM	0.495 mL	2.4751 mL	4.9502 mL
50 mM	0.099 mL	0.495 mL	0.990 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481