


Binankadsurin A

Chemical Properties

CAS No. : 77165-79-8
 Formula: C₂₂H₂₆O₇
 Molecular Weight: 402.44
 Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year
Actual storage temperature shall be subject to the COA.



Biological Description

| | |
|-------------|--|
| Description | Binankadsurin A is a lignan derived from Kadsura longipedunculata. |
|-------------|--|

Preparing Stock Solutions

| | 1mg | 5mg | 10mg |
|-------|-----------|------------|------------|
| 1 mM | 2.4848 mL | 12.4242 mL | 24.8484 mL |
| 5 mM | 0.497 mL | 2.4848 mL | 4.9697 mL |
| 10 mM | 0.2485 mL | 1.2424 mL | 2.4848 mL |
| 50 mM | 0.0497 mL | 0.2485 mL | 0.497 mL |

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481