


Tenacissimoside J

Chemical Properties

CAS No. : 1514895-60-3
 Formula: C₃₇H₅₈O₁₃
 Molecular Weight: 710.85
 Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year
 Actual storage temperature shall be subject to the COA.



Biological Description

Description	Tenacissimoside J is a steroidal glycoside that can be isolated from <i>Marsdenia tenacissima</i> .
-------------	---

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	1.4068 mL	7.0338 mL	14.0677 mL
5 mM	0.2814 mL	1.4068 mL	2.8135 mL
10 mM	0.1407 mL	0.7034 mL	1.4068 mL
50 mM	0.0281 mL	0.1407 mL	0.2814 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286 E_mail: info@targetmol.com Address: 34 Washington Street, Wellesley Hills, MA 02481