


(Z)-Tetradec-7-en-1-ol

Chemical Properties

CAS No. : 40642-43-1
 Formula: C₁₄H₂₈O
 Molecular Weight: 212.37
 Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year
 Actual storage temperature shall be subject to the COA.



Biological Description

Description	(Z)-Tetradec-7-en-1-ol acts as an activator for mitogen-activated protein kinase (MAP kinase).
Targets(IC50)	p38 MAPK

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	4.7088 mL	23.5438 mL	47.0876 mL
5 mM	0.9418 mL	4.7088 mL	9.4175 mL
10 mM	0.4709 mL	2.3544 mL	4.7088 mL
50 mM	0.0942 mL	0.4709 mL	0.9418 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481