

4-[3-(2,6-dihydroxyphenyl)-4-hydroxy-6-methyl-1-oxo-2,3-dihydro-1H-isoindol-2-yl]butanoic acid


Chemical Properties

CAS No. :

Formula: C₁₉H₁₉NO₆

Molecular Weight: 357.36

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year
Actual storage temperature shall be subject to the COA.



Biological Description

Description	4-[3-(2,6-dihydroxyphenyl)-4-hydroxy-6-methyl-1-oxo-2,3-dihydro-1H-isoindol-2-yl]butanoic acid is a high-purity biochemical reagent suitable for biochemical experiments and drug synthesis research.
-------------	---

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	2.7983 mL	13.9915 mL	27.983 mL
5 mM	0.5597 mL	2.7983 mL	5.5966 mL
10 mM	0.2798 mL	1.3991 mL	2.7983 mL
50 mM	0.056 mL	0.2798 mL	0.5597 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286 E_mail: info@targetmol.com Address: 34 Washington Street, Wellesley Hills, MA 02481