


22393-62-0_peak 2

Chemical Properties

CAS No. : 22393-63-1
 Formula: C₂₂H₁₉Br
 Molecular Weight: 363.29
 Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year
 Actual storage temperature shall be subject to the COA.



Biological Description

Description	22393-62-0_peak 2 exhibits some aspects of the antiestrogenic activity and other actions that may be connected to the estrogenic properties. A mixture of the (Z)- and (E)-isomers (Broparestrol, INN) is used in dermatology.
Targets(IC50)	Others

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	2.7526 mL	13.7631 mL	27.5262 mL
5 mM	0.5505 mL	2.7526 mL	5.5052 mL
10 mM	0.2753 mL	1.3763 mL	2.7526 mL
50 mM	0.0551 mL	0.2753 mL	0.5505 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481