


8-(1,1-Dimethylallyl)genistein

Chemical Properties

CAS No. :	651750-08-2
Formula:	C ₂₀ H ₁₈ O ₅
Molecular Weight:	338.359
Storage:	Powder: -20°C for 3 years In solvent: -80°C for 1 year Actual storage temperature shall be subject to the COA.



Biological Description

Description	8-(1,1-Dimethylallyl)genistein is a natural product for research related to life sciences. The catalog number is TN6335 and the CAS number is 651750-08-2.
-------------	--

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	2.9554 mL	14.7772 mL	29.5543 mL
5 mM	0.5911 mL	2.9554 mL	5.9109 mL
10 mM	0.2955 mL	1.4777 mL	2.9554 mL
50 mM	0.0591 mL	0.2955 mL	0.5911 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481