


Z 4212

Chemical Properties

CAS No. : 125483-43-4
 Formula: C₁₇H₂₄N₂O₆S
 Molecular Weight: 384.45
 Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year
Actual storage temperature shall be subject to the COA.



Biological Description

Description	Z 4212 is a biochemical.
Targets(IC50)	Others

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	2.6011 mL	13.0056 mL	26.0112 mL
5 mM	0.5202 mL	2.6011 mL	5.2022 mL
10 mM	0.2601 mL	1.3006 mL	2.6011 mL
50 mM	0.052 mL	0.2601 mL	0.5202 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Reference

Shao C, Tian Y, Dong Z, Gao J, Gao Y, Jia X, Guo G, Wen X, Jiang C, Zhang X. The Use of Principal Component Analysis in MALDI-TOF MS: a Powerful Tool for Establishing a Mini-optimized Proteomic Profile. *Am J Biomed Sci.* 2012;4(1):85-101. PubMed PMID: 22229059; PubMed Central PMCID: PMC3251008.

Tsuneki H, Ishizuka M, Terasawa M, Wu JB, Sasaoka T, Kimura I. Effect of green tea on blood glucose levels and serum proteomic patterns in diabetic (db/db) mice and on glucose metabolism in healthy humans. *BMC Pharmacol.* 2004 Aug 26;4:18. PubMed PMID: 15331020; PubMed Central PMCID: PMC517497.

Otero MJ, Iglesias T, Fuentes JA. Hypoalgesic action of bestatin analogues that inhibit central aminopeptidases, but not neutral endopeptidase. *Neuropeptides.* 1993 Sep;25(3):175-82. PubMed PMID: 8247255.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481