

Dihydromarein

Chemical Properties

CAS No. :

Formula: C₂₁H₂₄O₁₁

Molecular Weight: 452.4

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.

Biological Description

Description	Dihydromarein is available for research in the life sciences and its product number is TN6578.
-------------	--

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	2.2104 mL	11.0522 mL	22.1043 mL
5 mM	0.4421 mL	2.2104 mL	4.4209 mL
10 mM	0.221 mL	1.1052 mL	2.2104 mL
50 mM	0.0442 mL	0.221 mL	0.4421 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481