

## BMV Gag-(7-25)

## Chemical Properties

CAS No. : 312613-48-2


Formula: C96H170N44O24S

Molecular Weight: 2356.72

Keep away from moisture

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.



## Biological Description

Description	BMV Gag-(7-25) is a peptide enriched with arginine that possesses the ability to penetrate cells. It is utilized in research related to drug delivery and gene therapy.
Targets(IC50)	TMV

## Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	0.4243 mL	2.1216 mL	4.2432 mL
5 mM	0.0849 mL	0.4243 mL	0.8486 mL
10 mM	0.0424 mL	0.2122 mL	0.4243 mL
50 mM	0.0085 mL	0.0424 mL	0.0849 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286

E\_mail:info@targetmol.com

Address:34 Washington Street,Wellesley Hills,MA 02481