


(S)-PHA533533

## Chemical Properties

CAS No. : 437982-89-3  
 Formula: C<sub>19</sub>H<sub>22</sub>N<sub>4</sub>O<sub>2</sub>  
 Molecular Weight: 338.40  
 Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year  
 Actual storage temperature shall be subject to the COA.



## Biological Description

Description	(S)-PHA533533, a cyclin-dependent kinase 2 (CDK2) inhibitor, exhibits antitumor activity and has the capability to induce widespread expression of UBE3A in neurons of Angelman syndrome model mice.
Targets(IC50)	CDK

## Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	2.9551 mL	14.7754 mL	29.5508 mL
5 mM	0.591 mL	2.9551 mL	5.9102 mL
10 mM	0.2955 mL	1.4775 mL	2.9551 mL
50 mM	0.0591 mL	0.2955 mL	0.591 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286 E\_mail: info@targetmol.com Address: 34 Washington Street, Wellesley Hills, MA 02481