

Emerfetamab

Chemical Properties

CAS No. : 2250261-27-7

Formula:

Molecular Weight:

Storage: -20°C for 1 year
Actual storage temperature shall be subject to the COA.

Biological Description

Description	Emerfetamab (AMG 673) is a dual-specificity T-cell engager (bispecific T-cell engager) with an extended half-life that targets both CD33 and CD3. By redirecting T cells towards CD33+ cells, Emerfetamab induces T-cell mediated cytotoxicity against progenitor cells in acute myeloid leukemia (AML).
Targets(IC50)	Integrin

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481