

Trebananib

Chemical Properties

CAS No. : 894356-79-7

Formula:

Molecular Weight:

Storage: -20°C for 1 year
Actual storage temperature shall be subject to the COA.

Biological Description

Description	Trebananib (2xCon4C; AMG 386) is a monoclonal antibody targeting ANGPT2 (angiopoietin 2, Ang2).
Targets(IC50)	Tie-2

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481