

Emirodatamab

Chemical Properties

CAS No. : 2449199-61-3

Formula:

Molecular Weight:

Storage: -20°C for 1 year
Actual storage temperature shall be subject to the COA.

Biological Description

Description	Emirodatamab (AMG-427) is an anti-FLT3/CD3 bispecific T-cell engager with an extended half-life. It can induce high cytotoxicity in primary AML progenitor cells and enhance FLT-3 expression. Combining Emirodatamab with anti-PD-1 antibodies increases T-cell dependent cytotoxicity (TDCC). This compound is under investigation for use in relapsed or refractory AML.
Targets(IC50)	FLT

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only· Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481