


HIV-1 integrase inhibitor 11

Chemical Properties

CAS No. :	54030-51-2
Formula:	C ₈ H ₈ N ₄ O ₂ S
Molecular Weight:	208.24
Storage:	Powder: -20°C for 3 years In solvent: -80°C for 1 year Actual storage temperature shall be subject to the COA.



Biological Description

Description	HIV-1 integrase inhibitor 11 (compound 5) effectively inhibits HIV-1 integrase (HIV-1 integrase), with an IC ₅₀ of 125 μM.
Targets(IC ₅₀)	HIV Protease

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	4.8022 mL	24.0108 mL	48.0215 mL
5 mM	0.9604 mL	4.8022 mL	9.6043 mL
10 mM	0.4802 mL	2.4011 mL	4.8022 mL
50 mM	0.096 mL	0.4802 mL	0.9604 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481